




Measuring the LCLS few-fs xray pulse length

Andreas R. Maier^{1,3}


colaborators: Wolfram Helml^{2,3}, Wolfgang Schweinberger³, Justin Gagnon³, Vladislav Yakovlev³, Adrian L. Cavalieri⁴, Ivanka Grguras⁴, Stefan Düsterer⁵, Michael Meyer⁶, Paul Radcliffe⁶, Thomas Tschentscher⁶, Florian J. Gruener¹, Christopher Roedig⁷, Gilles Doumy⁸, John Bozek⁹, Ryan N. Coffee⁹, John Costello¹⁰ and Reinhard Kienberger^{2,3}

LMU¹ / TU² / MPQ³ / CFEL⁴/ DESY⁵ / XFEL⁶ / OSU⁷ / ANL⁸ / SLAC⁹ / DCU¹⁰


1. Measurement Method




2. Fighting the Jitter Problem




3. Single-Cylce THz Streaking





The Mission






GENESIS simulation of LCLS low-charge mode, 2fs FEL pulse (blue) and gaussian envelope (red). pulse close at saturation.

- LCLS is a pretty cool machine delivering tunable hard xray pulses
- we want to measure the actual pulse lenght
- for non-linear processes we need exact temporal pulse profile

required diagnostic:

- single-shot
- non-invasive
- ~fs resolution
- optical streak camera



Experimental Setup



CCD

- Experiment summer 2010 at LCLS AMO hutch
- high transmission of FEL pulse
- photoelectron pulse is a direct copy of FEL pulse profile

retardation plates

mirror

FEL

magnet


IR

LCLS FEL pulse

1keV - 2keV, few-fs pulse duration


IR laser

2 – 2.4 micron wavelength



Mapping time to energy






energy shift: $\Delta E \approx - p_e A(t)$ almost linear mapping from time to energy

2 micron IR laser, A=0.13 (red), GENESIS simulated 2keV 2fs FEL pulse (blue) $E_0 = 930 \text{ eV}$ initial photo electron kinetic energy, p_e : photo electron momentum, A(t): instantaneous vector potential


A. K. Kazansky, I. P. Sazhina, and N. M. Kabachnik: PRA 82, 033420 (2010) V. Yakovlev, J. Gagnon, N. Karpowicz, and F. Krausz: PRL 105, 073001, (2010)




Things are not that easy...



- summer 2010: only limited wavelength available: $2.4 \mu m 8 \text{ fs}$
- synchronization jiter FEL IR laser ~ 30-60 fs rms
- not single-shot analysis (yet)
- derive only typical pulse lengths







Analysis – The basic idea




4 fs FEL pulse







4 fs FEL pulse



4fs FEL pulse constructed from subspikes. central energy = 920 eV, λ = 2.4 μ m, A_{max} = .13 a.u. IR laser. 50 fs fwhm IR pulse length.







Analysis – The basic idea





4 fs FEL pulse






4 fs FEL pulse









Analysis – The basic idea



- correlate projected spectrum length and shift
- assumes pulse lengths of 2 fs and 10 fs




central energy = 920 eV / λ = 2.4 μ m / A_{max} = .13 a.u. / 20 fs rms jitter about the center of the 50 fs fwhm IR pulse





Somewhere around 4 fs fwhm








EXPERIMENTAL DATA analyze













Quantitative Analysis




- for every point in the correlation calculate all possible pulse lengths
- due to the unknown phase (jitter) pulse lengths are not unique





4 fs fwhm xray pulses






MEASURED:

Sort all possible pulse lengths calculate from measured spectra into histogram to get the typical pulse length.


This pulse length is an upper limit.

 4.4 ± 1.3 fs pulse width (fwhm)



Single-Spike, Single-fs Pulses






- Some spectra show single spikes.
- Corresponds to coherent single spike FEL pulses.

Back-of-the-envelope calculation:

- 4 fs mapped to 45 eV
- peak is 1 fs fwhm




Single-Cycle THz streaking



cross-correlation FEL with singlecycle THz pulse


Recent Experiment at FLASH (06/2011) by Adrian Cavalieri


- single-cycle THz pulse
- 300-500 fs rising edge
- shift of peak is time-of-arrival
- width of peak is pulse length
- limited by detector resolution
- independent of xray energy
- non-invasive




Summary











• 4 fs Pulses

- Shortest Xray pulses measured so far.
- Not yet single shot.



Single-shot measurement is possible.




Acknowledgment



Wolfi Helml and Paul Radcliffe for great help during data analysis.


Sven Reiche and Yuantao Ding for help with simulations and many discussions.

The whole LCLS team for support and that wonderful machine.




Summary





200 200 100 0 2 4 6 8 pulselength (fs fwhm)



Pulse length measurement with optical streak camera.

- 4 fs Pulses
- Shortest Xray pulses measured so far.
- Not yet single shot.

Single-shot measurement is possible.

Thanks for your attention.