


# Experiment Automation with a Robot Arm using the Liquids Reflectometer Instrument at the SNS

NEUTRON SCIENCES



Bogdan Vacaliuc, Piotr A. Zolnierczuk, Madhan Sundaram, Andre A. Parizzi, Candice E. Halbert, Michael C. Hoffmann, Gayle C. Greene, Jim F. Browning and John Francis Ankner




Neutron Sciences Directorate,  
Oak Ridge National Laboratory, Oak Ridge, Tennessee

## The SNS Liquids Reflectometer


The SNS Liquids Reflectometer is a horizontal-surface unpolarized instrument capable of operating in specular, off specular, and grazing-incident small angle neutron scattering (SANS) geometry.

### Operating parameters

- Wavelength range:  $2.5 \text{ \AA} < \lambda < 17.5 \text{ \AA}$
- Operating bandwidth:  $\Delta\lambda = 3.5 \text{ \AA}$
- 2D Position-sensitive detector resolution:  $1.3 \times 1.3 \text{ mm}^2$
- Solid-sample  $Q$  range:  $0.003 \text{ \AA}^{-1} < Q_z < 1.5 \text{ \AA}^{-1}$
- Air/liquid  $Q$  range:  $0.003 \text{ \AA}^{-1} < Q_z < 0.4 \text{ \AA}^{-1}$
- Minimum reflectivity (no hydrogen):  $5 \times 10^{-8}$
- Robotic sample changer with 18 wafer magazine



## Robot Control Network and System Design



The instrument control network provides access to motor controllers, sample environment control and detector data collection computers.

### Highlights

- Process variable distribution via NI DataSockets
- Shared memory system on control computer
- Kinematic mount enables sample holders to 'self-align'
- Fundamental 6-point robot motion for pick and put enables operation with sample magazine.



| Productivity Metrics                                                           |                 |                 |                                                                         |
|--------------------------------------------------------------------------------|-----------------|-----------------|-------------------------------------------------------------------------|
| Metric                                                                         | 2012-A          | 2013-A          | Notes                                                                   |
| Number of neutron runs that could have used the current robotic sample changer | 2696            | 3805            | Number of neutron runs that <i>did</i> use the robotic sample changer   |
| Number of days used to run experiments                                         | 30 days         | 51 days         |                                                                         |
| Average energy on target/experiment day                                        | 16.98 MWh/day   | 16.20 MWh/day   | <i>includes</i> downtime during experiments scheduled and non-scheduled |
| Average number of runs/energy on target                                        | 5.29 runs/MWh   | 4.61 runs/MWh   |                                                                         |
| Estimate average MWh used per sample                                           | 1.52 MWh/sample | 1.76 MWh/sample | based on 8 runs/sample (typical)                                        |
| Typical time of manual alignment scan                                          | 14 min/sample   | 32 min/sample   | Average time of automatic alignment                                     |
| Energy equivalent spent for alignment                                          | 165 kWh/sample  | 360 kWh/sample  |                                                                         |

| Future Projections                                                                   |                |  |
|--------------------------------------------------------------------------------------|----------------|--|
|                                                                                      | 2013-B         |  |
| Benchmarked "markers"-based alignment time                                           | 12 min/sample  |  |
| Energy equivalent estimate using "markers"-based alignment and accelerator at 875 kW | 135 kWh/sample |  |
| Expected number of runs/energy on target using "markers"-based alignments            | 5.29 runs/MWh  |  |

## Instrument Productivity Enhancement

### Robot Control

- Mechanical enhancements to ensure repeatability of table reassembly
- Vision system to verify position
- Satellite and control program
- PyDAS integration

| Instrument Control |                 |              |                                 |
|--------------------|-----------------|--------------|---------------------------------|
| Enable Robot       | Actual: DISABLE | Request:     | Status: <span>ATPOSITION</span> |
| Motor Positions    | UNDEFINED       |              | <span>OFFTARGET</span>          |
| Sample Select      | 4               | PICKUPPOS    | <span>ATPOSITION</span>         |
|                    |                 | RESTOREPOS   |                                 |
|                    |                 | MEMANDPICKUP |                                 |
|                    |                 | MEMPOS       |                                 |

### Automatic Alignment (solid samples)

- Uses neutron beam and data collection
- Aligns 3 axes of sample stage optimally
- Robust operation enables robot integration


## Lessons Learned

### Automation

- Initially, instrument operation can become "less efficient" as human-in-the-loop activity is replaced by algorithms
- Detailed timing analysis to identify bottlenecks and resolve them

### Vision System

- Consult with vendor FAE during setup
- Protect position-sensitive equipment from bumps



The authors wish to acknowledge valuable contributions made by Lloyd G. Clonts, Xiaodong Tao and Rick Reidel for their original work on the DAS control system. Robert Viola and "Cooper" at Square-One Systems for the design and programming of the robot arm. Brad Lokitz, Juan Pablo Hinestrosa, and Mike Kilbey for being the "first" users of the robotic experiment processing system.

Contacts: John Ankner ([anknerjf@ornl.gov](mailto:anknerjf@ornl.gov)), Jim Browning ([browningjf@ornl.gov](mailto:browningjf@ornl.gov))

Spallation Neutron Source at Oak Ridge National Laboratory, managed by UT-Batelle, LLC for the U.S. Department of Energy under Contract No. DE-AC05-00OR22725.

