



Micro-bunching measurement techniques in storage rings

Anke-Susanne Müller

INSTITUTE FOR BEAM PHYSICS AND TECHNOLOGY (IBPT)



Micro-bunching instability

- Complex and nonlinear dynamics in longitudinal phase space
- Can lead to bursts of very intense, coherent synchrotron radiation
- Observation of, e.g.,
 - longitudinal bunch profiles [“direct”]
 - emitted (coherent) radiation [“indirect”]
- Relevant timescales:
 - bunch length
 - dimension of micro-structures
 - distinguish between individual bunches
 - follow bursting behavior
 - observe slow changes (e.g., with current)
- Diagnostics requirements:
 - high resolution (ps) - high rate (500 MHz) - long term observation (secs - hrs)

Micro-bunching instability

- Complex and nonlinear dynamics in longitudinal phase space
- Can lead to bursts of very intense, coherent synchrotron radiation
- Observation of, e.g.,
 - longitudinal bunch profiles [“direct”]
 - emitted (coherent) radiation [“indirect”]
- Relevant timescales:
 - bunch length
 - dimension of micro-structures

(sub-) ps

 - distinguish between individual bunches
 - follow bursting behavior
 - observe slow changes (e.g., with current)
- Diagnostics requirements:
 - high resolution (ps) - high rate (500 MHz) - long term observation (secs - hrs)

Micro-bunching instability

- Complex and nonlinear dynamics in longitudinal phase space
- Can lead to bursts of very intense, coherent synchrotron radiation
- Observation of, e.g.,
 - longitudinal bunch profiles [“direct”]
 - emitted (coherent) radiation [“indirect”]
- Relevant timescales:

■ bunch length	(sub-) ps
■ dimension of micro-structures	
■ distinguish between individual bunches	500 MHz
■ follow bursting behavior	
■ observe slow changes (e.g., with current)	
- Diagnostics requirements:
 - high resolution (ps) - high rate (500 MHz) - long term observation (secs - hrs)

Micro-bunching instability

- Complex and nonlinear dynamics in longitudinal phase space
- Can lead to bursts of very intense, coherent synchrotron radiation
- Observation of, e.g.,
 - longitudinal bunch profiles [“direct”]
 - emitted (coherent) radiation [“indirect”]
- Relevant timescales:

■ bunch length	(sub-) ps
■ dimension of micro-structures	
■ distinguish between individual bunches	500 MHz
■ follow bursting behavior	10 ms

 - observe slow changes (e.g., with current)
- Diagnostics requirements:
 - high resolution (ps) - high rate (500 MHz) - long term observation (secs - hrs)

Micro-bunching instability


- Complex and nonlinear dynamics in longitudinal phase space
- Can lead to bursts of very intense, coherent synchrotron radiation
- Observation of, e.g.,
 - longitudinal bunch profiles [“direct”]
 - emitted (coherent) radiation [“indirect”]
- Relevant timescales:

■ bunch length	(sub-) ps
■ dimension of micro-structures	
■ distinguish between individual bunches	500 MHz
■ follow bursting behavior	10 ms
■ observe slow changes (e.g., with current)	sec - hrs
- Diagnostics requirements:
 - high resolution (ps) - high rate (500 MHz) - long term observation (secs - hrs)

Longitudinal phase space

■ Micro-bunching instability:



sim. charge distribution in long. phase space showing micro-structures





courtesy T. Boltz

sim. long. bunch profile
(very small variation!)

I.P.S. Martin et al., IPAC2015, MOPMA003




Indirect observation: emitted THz radiation




V. Judin et al., IPAC2010, WEPEA021

FFT of THz signal as a function of bunch current



V. Judin et al., IPAC 2014, MOPRO063



spectrum: dimension of typical structures

Classification of techniques - a try

THz techniques time domain

- spectrograms and information extracted from them
- “finger print” of instability

EO techniques far field

- detect photons from emitted THz pulse
- measure ‘wake fields’ (EOS) and bunch profiles (EOSD)

THz techniques frequency domain

- signature of sub-structures
- single-shot acquisition of radiation spectrum maps bunch dynamics

EO techniques near field

- directly detect bunch field
- measure wake fields (EOS) and bunch profiles (EOSD)

Classification of techniques - a try



THz techniques time domain

- spectrograms and information extracted from them
- “finger print” of instability

EO techniques far field

- detect photons from emitted THz pulse
- measure ‘wake fields’ (EOS) and bunch profiles (EOSD)

THz techniques frequency domain


- signature of sub-structures
- single-shot acquisition of radiation spectrum maps bunch dynamics

EO techniques near field

- directly detect bunch field
- measure wake fields (EOS) and bunch profiles (EOSD)

THz detectors - example ANKA

Quasi-optical broadband detector




room temperature
response time $< 200\text{ps}$
50 GHz up to 1 THz
based on schottky diode
ACST (acst.de)

Hot Electron Bolometer (NbN)



cryogenic (LHe)
response time $< 165\text{ps}$
200 GHz up to 4 THz
DLR high sensitivity [2]

YBCO detector



cryogenic (LN2)
response time $< 15\text{ps}$
30 GHz up to 2.5 THz
(KIT - IMS)
J. Raasch [1]

■ Resolve intensity of each bunch (minimal bunch spacing 2 ns)


[1] Thoma, P.; Raasch, J.; et al.; IEEE Trans. Appl. Supercond., vol.23, no.3, pp.2400206,2400206, June 2013

[2] Yu.P. Gousev, G.N. Gol'tsman, A.D. Semenov et al. Broadband ultrafast superconducting NbN detector for electromagnetic radiation. J. Appl. Phys., 75:3695, 1994.



courtesy M. Brosi

From THz signal to spectrogram

- Detector: must resolve bunches with a bunch spacing of 2 ns
 - e.g., zero-biased quasi-optical Schottky barrier diode



FFT




M. Brosi et al., IPAC 2016, TUPOR006

M. Brosi et al., arXiv:1605.00536


The faster way: “Snapshot” measurements with the KAPTURE system

■ KAPTURE: KARlsruhe Pulse Taking Ultra-Fast Readout Electronics

- Readout system compatible with: YBCO, HEB, and Schottky diode detectors
- Pulse amplitude and arrive time measurements with “mv” and “ps” accuracy
- Simultaneous acquisition of all buckets turn-by-turn in streaming mode
- Continuous acquisition for long observation time.
- Real-time data evaluation by GPUs




KAPTURE - Basic concept



courtesy M. Caselle




Frequency behaviour of CSR @ different bunch current



KAPTURE v2: compact+high sampling rate design


■ Target design:

- Very wide operating frequency from 0.2 to 3.5 GS/s per sampling channel (ANKA @ 1 GS/s)
- Modular design for multiple KAPTURE connections → to extend the sampling points over 4 samples
- PCB size reduced by factor 2, two KAPTURE connected to one readout card
- Form factor mechanically/electrically compatible with μ TCA system



courtesy M. Caselle

KAPTURE – 2. PCB technology




First prototype will be ready:
mid of September

courtesy M. Caselle


From THz signal to spectrogram with KAPTURE

- Cover beam current range with special filling pattern
- read out each bunch's time signal with KAPTURE




M. Brosi et al., IPAC 2016, TUPOR006

From THz signal to spectrogram with KAPTURE




- FFT for each bunch's signal:




M. Brosi et al., IPAC 2016, TUPOR006

From THz signal to spectrogram with KAPTURE



■ Sorted by bunch current:



M. Brosi et al., IPAC 2016, TUPOR006

Classification of techniques - a try



THz techniques time domain

- spectrograms and information extracted from them
- “finger print” of instability

EO techniques far field

- detect photons from emitted THz pulse
- measure ‘wake fields’ (EOS) and bunch profiles (EOSD)



THz techniques frequency domain

- signature of sub-structures
- single-shot acquisition of radiation spectrum maps bunch dynamics

EO techniques near field



- directly detect bunch field
- measure wake fields (EOS) and bunch profiles (EOSD)

Single-shot spectrometer

- Multi-band mm wave detectors
- Measurements at DIAMOND LS



A. Finn et al., IPAC 2016, MOPMR036

A. Finn et al., IPAC 2015, MOPTY080



Single-shot spectrometer

- Zooming into the bursts
 - Different behavior in different frequency ranges




courtesy J. Steinmann


Single-shot spectrometer

■ Integrated planar array of double-slot antenna- coupled YBCO detectors

3 mm by 3 mm four-element array




detector block



Poster J. Raasch et al. IBIC 2016, WEPG56

A. Schmid et al., ASC2016

single bunch in the Diamond LS filling pattern
simultaneously measured with all array channels




Classification of techniques - a try



THz techniques time domain

- spectrograms and information extracted from them
- “finger print” of instability



EO techniques far field

- detect photons from emitted THz pulse
- measure ‘wake fields’ (EOS) and bunch profiles (EOSD)



THz techniques frequency domain



- signature of sub-structures
- single-shot acquisition of radiation spectrum maps bunch dynamics

EO techniques near field


- directly detect bunch field
- measure wake fields (EOS) and bunch profiles (EOSD)

EO techniques at storage rings: first experiments


- SLS, FEMTO slicing experiment: EO sampling & EO spectral decoding measurements of CSR from 'sliced' bunch



- ANKA: asynchronous sampling



Phase sensitive detection of THz radiation with electro-optical femto-second sampling




S. Ibrahimkuty et al., J. Synchrotron Rad. (2011)18


A. Plech et al., PAC 2009 , TU5RFP026

EO in the far field: dynamics at large time scales

high-sensitivity photonic time-stretch
EO sampling



THz pulses shapes vs time



C. Evain, et al., arXiv:1607.08072 (2016)

C. Sz waj, et al., arXiv:1607.07168 (2016)

Classification of techniques - a try



THz techniques time domain

- spectrograms and information extracted from them
- “finger print” of instability



EO techniques far field

- detect photons from emitted THz pulse
- measure ‘wake fields’ (EOS) and bunch profiles (EOSD)



THz techniques frequency domain

- signature of sub-structures
- single-shot acquisition of radiation spectrum maps bunch dynamics




EO techniques near field

- directly detect bunch field
- measure wake fields (EOS) and bunch profiles (EOSD)


EO in the near field

Electro-Optical Spectral Decoding




EO in the near field

Electro-Optical Spectral Decoding



EO in the near field: bunch shape dynamics

- Single-shot EOSD measurements show dynamic sub-structures




SPONSORED BY THE



Federal Ministry
of Education
and Research


N. Hiller, et al., IBIC 2014, MOPD17

simulated phase space




courtesy J. Steinmann, P. Schönfeldt

KALYPSO: Karlsruhe Linear array detector for MHz-repetition rate Spectroscopy



■ Scientific goal: develop the “ideal” detector for EOSD measurements:

- High-repetition rate single-shot resolution
- Continuous data acquisition → turn-by-turn monitoring over $> 10^6$ turns
- Detect radiation in visible & near-infrared spectrum
- Real-time data analysis → FPGA / GPU heterogeneous DAQ



Poster M. Vogelsanget al.
IBIC 2016, WEPG07

■ Main technological challenges:

- Fast front-end electronics: GOTTHARD chip
- High density connections (wire-bondings, transmission lines, etc.)
- High-throughput DAQ system

courtesy L. Rota

KIT
Karlsruhe Institute of Technology

- 
- Poster L. Rota et al.
IBIC 2016, WEPG46**

2.7 MHz

GOTTHARD
"KIT-version"

ADC

InGaAs
sensor

DACs with
EEPROM

GOTTHARD
"KIT-version"

opening for
back illumination


- Aim: 2.7 MHz
(XFEL 5 MHz)
- Applications:
turn-by-turn long. &
transv. profiles

courtesy L. Rota

Single-shot bunch profiles with EO & KALYPSO

- Acquisition rate: 900 kHz - one image every 1 μs

courtesy P. Schönfeldt



time span (μs)

1E+01

1E+02

1E+03


1E+04

1E+05

Single-shot bunch profiles with EO & KALYPSO

- Acquisition rate: 900 kHz - one image every 1 μs

courtesy P. Schönfeldt



time span (μs)

1E+01

1E+02

1E+03


1E+04

1E+05

Single-shot bunch profiles with EO & KALYPSO

- Acquisition rate: 900 kHz - one image every 1 μs

courtesy P. Schönfeldt



time span (μs)

1E+01

1E+02

1E+03


1E+04

1E+05

Single-shot bunch profiles with EO & KALYPSO

- Acquisition rate: 900 kHz - one image every 1 μs

courtesy P. Schönfeldt



time span (μs)

1E+01

1E+02

1E+03


1E+04

1E+05

Single-shot bunch profiles with EO & KALYPSO

■ Acquisition rate: 900 kHz - one image every 1 μs

courtesy P. Schönfeldt



time span (μs)

1E+01

1E+02

1E+03


1E+04

1E+05

Single-shot bunch profiles with EO & KALYPSO

- Acquisition rate: 900 kHz - one image every 1 μs

courtesy P. Schönfeldt



time span (μs)

1E+01

1E+02


1E+03

1E+04

1E+05

Micro-bunching measurements with KALYPSO

Each vertical line corresponds to a single-shot profile measurement acquired at $f_{\text{rev}} = 2.7$ MHz




detail: development of sub-structures

courtesy L. Rota

Micro-bunching measurements with KALYPSO

■ Synchronous measurements with KALYPSO and KAPTURE



courtesy L. Rota

Summary

- Micro-bunching measurement techniques allow to study the complex and nonlinear dynamics in longitudinal phase space
- Observation of, e.g.,
 - longitudinal bunch profiles [“direct”]
 - emitted (coherent) radiation [“indirect”]
- Diagnostics requirements:
 - high resolution (ps) - high rate (500 MHz) - long term observation (secs - hrs)
- Recent developments based on:

**THz techniques
time domain**

**EO techniques
far field**

**THz techniques
frequency domain**

**EO techniques
near field**

- Next steps?
 - combination of methods?

Acknowledgements

- T. Boltz, M. Brosi, F. Rämisch, M. Schuh, P. Schönfeldt, M. Schwarz, J. Steinmann, and the rest of the ANKA THz team
 - J. Raasch, A. Schmid (KIT-IMS), M. Caselle and L. Rota (KIT-IPE)
 - S. Bielawski (Université de Lille) and N. Hiller (PSI)
- ... for instructive discussions and for providing material and plots

